

BRONCHOSCOPE OP-40 / OP-50

GREAT VERSATILITY

- High resolution
- Easy handling
- Working channel and suction



OPTOMIC OP-40 and OP-50 flexible endoscopes are especially designed to facilitate passage through the trachea to permit diagnostic and therapeutic procedures.

BRONCHOSCOPE **OP-40 / OP-50**



FIBROLUX 250
LIGHT SOURCE



FIBROLUX 150
LIGHT SOURCE



FIBER OPTIC CABLE



INSTRUMENTS



LEAKAGE TEST DEVICE



WORKING CHANNEL

	OP-40	OP-50
Optical system		
Field of view	95° forward	100° forward
Depth of field	5-50mm	3-50mm
Minimum visible distance	5mm	3mm
Insertion tube		
Outer diameter	3.9mm	4.9mm
Distal diameter	3.9mm	4.9mm
Working length	600mm	600mm
Total length	850mm	850mm
Range of tip bending		
Up	160°	160°
Down	100°	100°
Biopsy channel		
Working channel inner diameter	1.4mm	2.0mm

OPTOMIC ESPAÑA S.A.

Madroño, 4. E-28770 Colmenar Viejo, Madrid, Spain
Tel: +34 902 22 11 77 / +34 91 803 57 14 • Fax: +34 902 22 11 78
export@optomic.com • www.optomic.com



Specifications are subject to change without notice

