



Ward Trolleys & Carts

www.francehospital.com

francehospital®
trolleys for healthcare









We begin each product design and every customer solution with clear criteria born of sound research into client needs. Our research summaries gather and interpret studies of workplace issues, offering a wide range of solutions for any specific applications. Constant communication with our clients allows us to collect their requests and to improve constantly our products which are therefore made according to healthcare's high standards of hygiene, safety and durability.

Ward Trolleys & Carts



Index

1.	<i>PRECISO</i>	6
2.	<i>PERSOLIFE</i>	12
3.	<i>PERSOCAR</i>	16
4.	<i>PERSOCARCOLOR</i>	20
5.	<i>COLORLIGHT</i>	22
6.	<i>SHANGO</i>	24
7.	<i>PERSOCAR LIGHT</i>	26
8.	<i>ERGODOC</i>	28
9.	<i>Utility Carts</i>	30
10.	<i>Internal Organzier for Drawers</i>	34
11.	<i>Accessories for all the trolleys</i>	36

<i>Cart selection/Premium solution</i>	<i>Top view</i>	<i>Body sizes</i>	<i>Anaesthesia</i>	<i>Emergency</i>	<i>Therapy</i>	<i>Dressing</i>
Preciso		830x649 mm	 PRE9AV-G-ANE p. 7		 PRE9XAPS-G-FAR p. 11	 PRE8AER-G-CHI p. 7
Persolife		660x550 mm		 PCL p. 12		
Persocar		955x550 mm 660x550 mm	 PCMA2 p. 18			 PCM2 p. 18
Persocarcolor		660x550 mm	 PCM2-V-S p. 21			 PCM2-V p. 21
Colorlight		830x545 mm	 CLM5CP-I p. 23			 CLM4CP p. 22
Shango		660x510 mm				 SHM p. 25
Persocar Light		455x530 mm	 PCMA3 p. 26	 PCL3C p. 27		 PCM3 p. 26
Ergodoc		920x650 mm				

<i>Visit</i>	<i>Suture</i>	<i>Treatment</i>	<i>Logistics</i>	<i>Paediatric/Broselow</i>	<i>Technology</i>
 <i>PRE9XAB-S-VIS</i> p. 7		 <i>PRE9AA-G-TER</i> p. 6	 <i>PRE9X-AAA-S</i> p. 9	 <i>PRE-GP</i> p. 7	 <i>PRE9X-APR-SPC-RIA</i> p. 7
				 <i>PCL-V-B</i> p. 13	
	 <i>PCM42</i>	 <i>PC2/PCM4</i> p. 16			 <i>PCM2-S</i> p. 18
	 <i>PCM22-V</i>	 <i>PC2-V</i> p. 20			 <i>PCM2-V-S</i> p. 21
	 <i>CLM-SU</i>	 <i>CL5CP</i> p. 23		 <i>CL-GP</i> p. 23	
		 <i>SH4G-V</i> p. 25		 <i>SH-S</i> p. 25	 <i>SH4G-G</i> p. 24
		 <i>PC3</i> p. 26		 <i>PC35</i> p. 27	
 <i>ED1</i> p. 29					

PRECISO

Modular system for a wide range of ward trolleys



PRE9AAA-G-TER



PRE9AAA-G-MED

When creativity meets innovation the result is amazing!

Assembled structure with angulars and lateral walls in extruded anodized aluminium, back wall in multilayer aluminium-polypropylene-aluminium.

Spacious worktop (DIN A2 sheet) in Polycarbonate/ABS injection moulded plastic with the innovative and unique handles on 3 sides.

Base in Polycarbonate/ABS injection moulded plastic with all-around bumpers.

Lateral extendable writing surface

Twin wheels Ø 125 mm, 2 with brake

On request: lateral additional drawer with independent key-lock, usable for laptop PC or narcotics.



Colour choice - neutral, pastel and electric colours:



PRExAEV...



PRExAEG...



PRExAER...



PRExAEB...



PRExAPS...



PRExAPR...



PRExAPB...



PRExAPV...



PRExAAA...

PRECISO with drawers



Unique drawers in Polycarbonate/ABS injection moulded plastic with frontal handle complete with coloured ID tag and label holder.

Dim.: ext. 610x470 mm

Int.: 540x410 mm

Available in 4 heights:

- 76 mm (3") - 1 module
- 156 mm (6") - 2 modules
- 236 mm (9") - 3 modules
- 316 mm (12") - 4 modules

The drawers can be used on ball-bearing rails - full-extending, self-closing and with anti-tilting feature - on a frame.

Load capacity 45 kg each drawer

The drawer can be easily removed to exchange and cleaning





As an alternative the drawers can be used on the plastic ISO-rails, like the standard ISO-Modules (baskets and trays).

Load capacity 25 kg each drawer

You can easily change the drawer configuration.



ID tag available in 9 colours

Ruby Red	 Mod. IDR	Cadmium Yellow	 Mod. IDG	Lily Violet	 Mod. IDVI	Traffic Orange	 Mod. IDAR	Pale Pink	 Mod. IDRO
Mint Green	 Mod. IDV	Gentian Blue	 Mod. IDB	Vaio Gray	 Mod. IDGR	Cream White	 Mod. IDBI		

PRECISO with ISO-rails and ISO-Modules

Closing by rolling shutter or doors



Mod. CP-1M / CP-2M / CP-3M

*Drawers 76 mm - 156 mm - 236 mm
(3-6-9 ISO-rails)*



Mod. VP050/100/200 - CP100/200

*ISO-trays and -baskets
600x400x h 50 mm - 100 mm - 200 mm
(3-5-8 ISO-rails)*



Mod. PRE-SV

*Vertical dividers for patient files
(min. 13 ISO-rails for headroom
h. 30 cm DIN A4)*



Mod. CA-3M / CA-4M / CA-5M

*Bin-Cassettes for Medicine-Bins
(9-12-15)*



Mod. PRE-R

*Internal height adjustable shelves
(1 ISO-rail)*




Mod. VP-CM & UD-CM

*OTTAVIO: 8 Medicine dispensers on
a tray for daily therapy
(3 ISO-rails)*



Mod. PRE-TEC

*Framed patient file holder for patient
charts A4 dim. 330x250 mm or X-ray
dim. 465x390 mm
(min. 12 ISO-rails for DIN A4 charts)*



Mod. VA-BIN & BIN

*Bin-Trays for Medicine-Bins
(4 to 6 for each tray)
(3-6 ISO-rails)*



PRE9APV-A



PRE9XAAA-S



PRE9XAAA-S

**>>> see ISO-rails capacity on the
Technical Specifications table**

Multiple accessories and different configurations to meet every need

*** For accessories see page 34-39



PRE9XAPR-SPC-RIA



PRE9AAA-G-INF



PRE9XAPS-G-DRE

Technical specifications							
Mod.	N°7	N°8	N°8X	N°9	N°9X	N°10	N°10X
Product codes							
- Drawer on rails	PRE7xxx-G	PRE8xxx-G	PRE8Xxxx-G	PRE9xxx-G	PRE9Xxxx-G	PRE10xxx-G	PRE10Xxxx-G
- Drawer on ISO	PRE7xxx-I	PRE8xxx-I	PRE8Xxxx-I	PRE9xxx-I	PRE9Xxxx-I	PRE10xxx-I	PRE10Xxxx-I
Drawer headroom	7 modules	8 modules	8 modules	9 modules	9 modules	10 modules	10 modules
ISO-rails	21	24		27		30	
Dimensions (mm)	830x649						
Worktop h. (mm)							
- drawers	840	920	1020	1000	1100	1080	1180
- doors		940		1020		1100	
Wheels	Twin 4x Ø 125 mm – 2 with brake						
Options on request	Doors	Doors	Rolling shutter / lateral drawer	Doors	Rolling shutter / lateral drawer	Doors	Rolling shutter / lateral drawer



PRE9XAPS-G-FAR



PRE9AAA-G-BRO



PRE9APB-G-GEN



PRE9AAA-S-LOG

PERSOLIFE Emergency trolley

Persolife is the most compact emergency trolley on the market. Thanks to its ergonomic and smart design, it stores all the necessary equipment and devices, allowing the contemporaneous access to three operators.

Available in two versions, with structure and drawers:

1. in powder coated steel with anti-microbe and anti-bacteria treatment, in dark-grey
2. in stainless steel AISI 304

Quickness and efficiency

- Drawers in several sizes, mounted on fully-extending top-quality rails
- Multiple drawer organizers and dividers for easy, quick access to drugs and supplies
- Colored handles with label holder, for quick identification of content
- Breakable, numbered seal ensures safety of drawer content



PCL



PCL-V-KL



PCLM



PCL-V-K

Broselow colour coded cart, for paediatric emergency





PCL front

Specific incorporated compartments to reduce footprint

- Suction pump holder
- Tilt-out bins
- Cardiac board
- Extendable writing surface
- Closing by numbered breakable plastic seal



PCL back

Everything in its place

- IV pole in stainless steel with two hooks, height adjustable
- Defibrillator plate, 360° swivel with locking mechanism, stainless steel arm and size adjustable Plexiglas plate with fixing strips.
- Oxygen bottle compartment, from 3 to 7 liters
- Storage space for catheters / probes
- Waste container with removable internal bin



When cardiac board is not required, or preferred on a lateral side:
closing with knob and plastic seal.

*** For accessories see page 34-39

Technical specifications						
Product codes	PCL	PCL-V	PCL-V-K	PCL-V-KL	PCL-V-B	PCLM
	In stainless steel	In powder coated steel	Without cardiac board	With lateral cardiac board	With 8 Broselow code drawers	With overbridge
Drawer configurations	2x h. 75 mm – 2x h. 155 mm – 1x h. 235 mm					3x h. 155 mm – 1x h. 235 mm
Dimensions (mm)	780x690x1250 H mm					H 1600
Body / Worktop h.	600x550 / 1020 h					
Weight (kg)	63	65	59	62	70	68
Wheels	Ø 150 mm, 2 antistatic, 2 with brake – double ball bearing					

- Drawers available in 4 heights:
 - h 75 mm (int. 360x470x45 h mm) – 1 module
 - h 155 mm (int. 360x470x135 h mm) – 2 modules
 - h 235 mm (int. 360x470x205 h mm) – 3 modules
 - h 315 mm (int. 360x470x290 h mm) – 4 modules
- Centralized closing by numeric breakable plastic seal (on request by key lock and seal)
- Coloured strips for content drawers identification
- Spacious worktop in moulded plastic with rim on 3 sides in light-grey colour
- Built-in push handle



PCL-V

PERSOCAR

Multifunction trolley in stainless steel AISI 304

The original design and the smart constructive solutions make the success of the Persocar line. It has structure and drawers in stainless steel AISI 304 and moulded plastic mint-green sides; this combination guarantees safety, reliability and durable life cycle. Extremely flexible and multifunction, it can be configured for specific uses: medication, dressing, anaesthesia, isolation, suture, and many others.

Persocar line is available in two dimensions:

- > 600x550 mm (mod. PC2)

- > 950x550 mm (mod. PC4)

to meet every needs of space.



PC2



PC4 with accessories



- Drawers available in 4 heights:
 - h 75 mm (int. 360x470x45 h - 720x470x45 h) – 1 module
 - h 155 mm (int. 360x470x135 h - 720x470x135 h) – 2 modules
 - h 235 mm (int. 360x470x205 h - 720x470x205 h) – 3 modules
 - h 315 mm (int. 360x470x290 h - 720x470x290 h) – 4 modules
- Stainless steel rails with full extraction
- Centralized closing by key lock
- Spacious worktop in moulded plastic with rim on 3 sides in light-grey colour
- Built-in push handle
- Lateral extendable work surface
- 3 tilt-out bins in see-through plastic
- MAXIWALLY knee-opening waste bin, with internal extractable bucket 20 l.



PC2



PC4-APF

Multiple accessories and different configurations to meet every need



PCMA2



PCM2



PCM23

Technical specifications

Product codes	PC2 / PCM2	PC4 / PCM4
	Drawers / With overbridge	Drawers / With overbridge
Drawer headroom	9 modules	9 modules
Dimensions (mm)	850x620x1070 / 1700 with overbridge	1180x620x1070 / 1700 with overbridge
Body / Worktop h.	600x550 / 1020 h	955x550 / 1020 h
Weight (kg)	49 (basic configuration)	71 (basic configuration)
Wheels	4 Ø 150 mm, 2 with brake	



PC2: built-in catheter holder on back side

Dozens of accessories to configure your trolley

*** For accessories see page 34-39



PCM4-S



PC4D

PERSOCARCOLOR

Multifunction trolley in powder coated steel



The **Persocarcolor** line has structure in powder coated steel with anti-microbe and anti-bacteria treatment, and moulded plastic mint-green sides. Extremely flexible and multifunction can be configured for specific use: medication, dressing, anaesthesia, isolation, and many other.

- Drawers available in 4 heights:
 - h 75 mm (int. 360x470x45 h mm) – 1 module
 - h 155 mm (int. 360x470x135 h mm) – 2 modules
 - h 235 mm (int. 360x470x205 h mm) – 3 modules
 - h 315 mm (int. 360x470x290 h mm) – 4 modules
- Centralized closing by key lock or
- Worktop in moulded plastic with rims on 3 sides
- Built-in push handle
- Lateral extendable work surface
- 3 tilt-out bins in see-through plastic
- MAXIWALLY knee-opening waste bin 20 l.



Technical specifications

Product codes	PC2-V	PCM2-V	PCM21-V	PCM22-V
	4 drawers	Overbridge & dispenser	Overbridge & trays	Overbridge & suture
Drawer headroom	9 modules			
Dimensions (mm)	850x620x1070 / 1700 with overbridge			
Body	600x550 / 1020 h			
Weight (kg)	49 (basic configuration)			
Wheels	4 Ø 150 mm, 2 with brake			

Built-in catheter holder on back side.

*** For accessories see page 34-39

COLORLIGHT

Multifunction trolleys in opal-blue powder coated steel

Technical specifications		
Mod.	With drawers	With overbridge
Product codes	CLxCP	CLMxCP
Drawer headroom	9 modules	9 modules
Dimensions (mm)	1040x555	
Body / Worktop h.	950	
Weight (kg)	70 (basic)	
Wheels	Twin 4x Ø 125 mm – 2 with brake	

The **Colorlight** line is available with unique drawers in Polycarbonate/ABS injection moulded plastic running on ball-bearing rails, full-extending, self-closing and with anti-tilting feature.

- Epoxy powder coated monoframe structure, with anti-microbe and anti-bacteria treatment (resistant for the entire life cycle of the trolley)
- Spacious worktop in moulded plastic with rim on 3 sides in light-grey colour
- Built-in push handles on both sides
- Lateral extendable working surface
- 3 tilt-out bins in see-through plastic
- MAXIWALLY knee-opening waste bin, with internal extractable bucket 20 l.



CLM4CP

Drawers

The drawers are complete with frontal handle complete with coloured ID tag and label holder.

The drawer can be easily removed to exchange and cleaning.

Load capacity 45 kg each drawer

Dim.: ext. 610x470 mm – Int.: 540x410 mm.

Available in 4 heights:

- 76 mm (3")
- 156 mm (6")
- 236 mm (9")
- 316 mm (12")

>>> see page 34-39 for the drawer internal organizers.



CLM5CP-I



CL5CP



CL4CP



SHANGO

Compact, with laptop computer drawer, in powder coated steel



SH4C-G

Shango combines the features of a traditional ward trolley with the innovative daily use of a laptop computer to input patient data, check prescription, visualizing X-rays, and more. The lateral laptop computer drawer has space for cables and battery and independent key lock.

With its pleasant design and compact size, **Shango** is extremely versatile.

- Powder coated steel structure, in light-grey colour, with anti-bacteria and anti-microbe treatment
- Powder coated steel drawers available in 4 colours:
 - Dark-grey
 - Melon-yellow
 - Mint-green
 - Sky-blue
- Spacious moulded worktop, in light-grey plastic, with rims on 3 sides
- Built-in push handle
- Drawers in four different heights and centralized key lock (on request: anti tilt system and electronic digit pad lock):
 - h 75 mm (int. 360x465x70 h mm) – 1 module
 - h 155 mm (int. 360x465x150 h mm) – 2 modules
 - h 235 mm (int. 360x465x230 h mm) – 3 modules
 - h 315 mm (int. 360x465x310 h mm) supplied with rails for hanging files - 4 modules
- 3 tilt-out bins in see-through plastic



SHM



SH-S



SH4G-V



SH-B-S

Dozens of accessories to configure your trolley

With overbridge and accessories rails height adjustable, various waste bins, IV pole, syringes distributor, stainless steel shelves, moulded plastic trays, gloves box holder, kidney dish holder, disinfectant instrument tray, laptop computer holder, catheter holder, and many others.

*** For accessories see page 34-39



Technical specifications				
Product codes	SH-G	SH	SH-V	SH-B
Drawers colours	Dark-grey	Melon-yellow	Mint-green	Sky-blue
Drawer headroom	9 modules			
Dimensions (mm)	660x510x1100 / 1730 (with overbridge)			
Body / Worktop h.	600x510 / 1050 h			
Weight (kg)	62 (basic)			
Wheels	Twin 4 Ø 125 mm, 2 with brake			

PERSOCAR LIGHT

Multifunction light trolley in stainless steel

Extremely compact and multifunction, can be used everywhere and configured for several purposes: emergency, medication, dressing, anaesthesia, isolation, suture, and many others.

- Drawers available in 4 heights:
 - h 75 mm (int. 360x470x45 h) – 1 module
 - h 155 mm (int. 360x470x135 h) – 2 modules
 - h 235 mm (int. 360x470x205 h) – 3 modules
 - h 315 mm (int. 360x470x290 h) – 4 modules
- Coloured strips for drawers content identification
- Worktop in moulded plastic with rims on 3 sides in light-grey colour
- WALLY knee-opening waste bin, with internal extractable bucket 8 l.

*** For accessories see page 34-39



PC3



PCM3

Technical specifications		
Product codes	PC3 / PCM3 / PCL3	PC35
	Drawers / With overbridge / Emergency	Drawers
Drawer headroom	9 modules	5 modules
Dimensions (mm)	685x600x1000 /1630 with overbridge	550x590x980
Body / Worktop h.	455x530 / 970 h	455x530x970
Weight (kg)	48 (basic configuration)	25
Wheels	4 Ø 125 mm, 2 with brake	



PC35



PCL3C front

For the emergency

In addition on the **Persolife Light**:

- Centralized closing system by numeric breakable plastic seal
- Oxygen bottle holder (from 3 and 5 litres)
- Catheter holder
- Cardiac board
- Defibrillator platform, size adjustable, 360° routable
- IV pole in stainless steel with 2 hooks, height adjustable
- Antistatic wheels

See also the **PERSOLIFE** emergency trolley page 12



PCL3C back

ERGODOC Patient file trolley



ED1-back

With the ward round begin one of the most important moment in the daily routine of a ward.

Ergodoc provides a modern and innovative solution for this need, providing greater productivity (storage space for gloves, disinfectant and a space for any other items) and more efficiency (ergonomic design). Moreover, it is easy to manage (spacious and extremely lightweight).





ED1-front



Accessories on request:

hanging files: DIN A4 (330x250 mm) – C1-FH,
Foolscap, X-ray size dim. 465x390 – C3-FH

- Structure in epoxy powder coated steel – fronts in vacuum moulded plastic
- Spacious container 477x477 mm on top, equipped with a frame to hold hanging files from A4 to foolscap and large X-Ray size. Capacity, about 35 files. Closed by a working surface with key-lock.
- Purple worktop (580x480 mm), in polyurethane. When opened it is a large working surface, with a perfect ergonomic height of 105 cm.
- Integrated handles on both sides
- Two lateral openings 18x16 cm

- Bottle holder for disinfectant with a fixation strap
- Spacious tilt-out bin
- Big lateral opening for various items 38x32x9 cm
- Specific opening to contain 3 boxes of gloves (S-M-L)
- Frontal large drawer holding up to 15 files A4 plus space on the side



Technical specifications	
Product Codes	ED1
	With drawer
Dimensions (mm)	920x650x1040 h
Weight (kg)	43
Wheels	4 Ø 125 mm, 2 with brake

Stainless steel utility carts

Entirely in stainless steel AISI304, these multifunction trolleys are from a functional and hygienic point of view, the ideal solution for daily transport inside the wards.



CS962



CS963



CS962-1C
1 drawer h 155

With moulded shelves

Technical specifications		
Products codes	CS962	CS963
Shelves	2 shelves	3 shelves
Dimensions (mm)	1000 x 690 x 1005	
Materials	stainless steel AISI 304	
Wheels	4 Ø 125 mm, 2 with brake	



CS962-2CP
2 drawers h 155



CS962-2C
2 drawers h 155 - 235



CS643RA

With removable shelves

Technical specifications			
Products codes	C3P000-06-4	C3P000-08-5	C3P000-10-5
Shelves	3 (available also with 2)		
Dimensions (mm)	680 x 490 x 950	880 x 590 x 950	1080 x 690 x 950
Materials	stainless steel AISI 304		
Wheels	4 Ø 125 mm, 2 with brake		



UCM44

Specific solutions

MEDONE multipurpose trolley two shelves (400x400) with railing to use in a various situations: medication, therapy, personal care of the patient, etc. We can produce custom made models to suit all needs. On request with antistatic wheels.

Tailor-made solution on demand

Technical specifications	
Products codes	UCM44
Dimensions (mm)	500x470x1010 h
Materials	stainless steel AISI 304
Wheels	4 Ø 125 mm, 2 with brake

PERMODUL Trolleys and shelving units



PERMODUL system combines the lightness of the aluminium posts with the strong polymer of the shelves. The two elements can easily be assembled without tools thanks to the ingenious fastening system. The specific bases allow to compose both trolleys and fixed shelving units. The original ergonomic push handle protects the hands and helps driving the trolley.

Several accessories complete the multipurpose trolleys.




PMC52



PMC53

Technical specifications	
Product codes	PMC52 & PMC53
Dimensions (mm)	920x500x1040 h
Capacity (kg)	160
Materials Shelves Posts	polymer anodised extruded aluminium
Wheels	4 Ø 125 mm, 2 with brake



MOSYS Wire chromed utility carts

The MOSYS carts can be used in several environments: storage, exhibition areas, archives, shops, cellars, cold rooms, washing room, wards - for all transport and logistical needs.

Available with shelves in multiple dimensions: from 400x600 mm to 600x1800 mm

*The MOSYS system is easy-to-assemble, without tools.
The shelves are height adjustable each 100 mm*



UC1




UC2



UC3

MOSYS accessories



dividers - D



edge fences - L



plastic bins/drawers - BOX




label holders - BM




Internal organizers for drawers

TRAYS

to optimize the drawer organization




On demand: removable and adjustable dividers and unique label holders:




The mauve bins are available in several sizes to mix depending on the needs:

Length 400 mm – Width 82 mm	
H. 58 mm	H. 138 mm
VSC6-1M	VSC6-2M
	
Length 400 mm – Width 125 mm	
H. 58 mm	H. 138 mm
VSC4-1M	VSC4-2M
	
Length 260 mm – Width 125 mm	
H. 58 mm	H. 138 mm
VSC42-1M	VSC42-2M
	



SOPPY

Comb-shape dividers, in transparent methacrylate, to combine following the specific requirements. They are suitable for all the trolleys.



PCSOP

For the emergency trolley PERSOLIFE and PERSOCAR trolleys – specific moulded trays complete with removable dividers:

PCSOP5 for ampoules		PCSOP48I for 4 rows and 16 dividers and 8 moulded spaces		PCSOP2I 2 rows and 6 dividers	
PCSOP3I 3 rows and 11 dividers		PCSOP4I 4 rows and 16 dividers		PCSOP2IA 2 rows and 6 dividers -height 130 mm	

Accessories for all the trolleys



Back accessory rail



I.V. pole



Stainless steel catheter basket



Basket in stainless steel wire mesh



Light bag-holder - 10 l.




Lateral accessory rail



Twin Trash bins with support



MAXIWALLY trash bin - large size



WALLY trash bin - small size



Trash bag holder




Mod. PCMSC

2 plastic bins with lids



Mod. PCMBC54

Syringe/needles dispenser



Mod. PCMMI2

Stainless steel tray




Mod. VA

Plastic tray – large size



Mod. VA-D

Plastic Tray



Mod. BIN42-1M/2M

Supply bins



Mod. SPG-T2/T3

Glove box holder (for 1 or 3 glove boxes)



Mod. PCMF

I.V. hook



Mod. PCMA

Kidney dish holder



Mod. PCMVD

Instrument disinfectant basin

Accessories for all the trolleys



Disinfectant dispenser




Cups dispenser



Sharp container holder



Instrument hook



Tape dispenser



Laptop PC holder



Stainless steel basin Ø 150x70 h mm



Stainless steel bin 280x200x550 h mm



Night light with clamp




Vertical accessory tube Ø 25x300 h mm



Mod. TPC-CL

Plastic box for 10 hanging charts A4 size



Mod. PCxSAC

Rails pair for hanging charts A4 size -
drawer h. 315 mm



Mod. PCVSPCAN

Steel frame for hanging charts A4 size



Mod. TPD

Steel frame for X-rays



Mod. C1-FH & C3-FH

Hanging charts A4 & X-rays/foolscap sizes



Mod. DW-SP

Double hook for Planetta



Mod. PI-xxx

Stainless steel plate to protect worktop



Mod. UCPF

Bottle holder - 4 holes Ø 80 mm



Mod. CL1000SG

Electronic keypad lock



For paediatric wards: fantasy graphics



FRANCEHOPITAL S.A.S.

*Z.I. Ouest Rue Georges Besse
B.P. 50030
F-67151 ERSTEIN Cedex / FRANCE
Tél. +33 03 88 59 87 87
Fax +33 03 88 98 04 44
francehopital@francehopital.fr*

*Export Office
via Bolzano, 2
I-39011 Lana / BOLZANO - ITALY
Tel. +39 0473 55 26 11
Fax +39 0473 55 26 99
italia@francehopital.com
export@francehopital.com*